# IR SYSTEM

## AM0180 - AM0500

### INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 12/2021

Revision: B

EC REP

**EC Representative**  
Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

**Manufacturer**  
Adept Medical Ltd  
2-6 McDonald St,  
Morningside  
Auckland 1025, New Zealand

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



### AM0180 (STARBoard + Extension Tray)

	<b>A</b>	STARBoard x 1
	<b>B</b>	Extension Tray x 1
	<b>C</b>	Arm Support Pad x 1

### AM0500 (IR Shield)

	<b>D</b>	X-Ray Shield x 1
	<b>E</b>	Daggerboard x 1

## Setup

(Instalación / Installation / Configurazione / Einrichtung / Installatie)

### STARBoard + Extension Tray

Check friction joints and adjust if necessary



- / Compruebe las juntas de fricción y ajústelas si es necesario
- / Contrôlez les articulations à friction et réglez-les si nécessaire
- / Controllare l'attrito delle articolazioni e, se necessario, regolare
- / Reibverbindungen prüfen und ggf. einstellen
- / Controleer frictieverbindingen en stel ze indien nodig af

Tighten  Loosen

- / Apretar
- / Serrer
- / Serraggio
- / Festziehen
- / Vastdraaien


 Aflojar

- / Desserrer
- / Allentamento
- / Lösen
- / Losdraaien


Refer to page 9

Consulte la página 9 / Reportez-vous à la page 9 / Vedere pag. 9 /  
Siehe Seite 9 / Zie pagina 9


3


**Taping (Optional)**

Encintado (opcional) / Contention par bande adhésive (facultatif) / Fissaggio con nastro (opzionale) / Fixieren (Optional) / Vastplakken (optioneel)


**Hand Positioning Options**

Opciones de posicionamiento de las manos / Options de positionnement de la main / Opzioni di posizionamento delle mani / Optionen für die Handpositionierung / Opties voor handpositie


**IR Shield (AM0500)**


1

2

3

4


## Spare Parts

(Piezas / Pièces de rechange / Parti di ricambio / Ersatzteile / Reserveonderdelen)

### Product Code

(Código de producto / Code produit / Codice prodotto / Produktcode / Productcode)

	STARBoard + Extension Tray	AM0180
	IR Shield	AM0500
	Arm Support Pad	AM0171

### Component Codes

(Códigos de componentes / Codes des composants / Codici dei componenti / Komponentencodes / Onderdeelcodes)


	Extension Tray	M3334
	X-Ray Shield	AM0630
	Daggerboard	M3677

### Maintenance instructions:


Instrucciones de mantenimiento: / Consignes d'entretien: / Istruzioni per la manutenzione: / Wartungshinweise: / Onderhoudsinstructies:

adeptmedicaltraining.com  
/product-maintenance-starsystem


## STARBoard

### STARBoard Component Codes

Códigos de componentes de STARBoard / Codes des composants du STARBoard / Codici dei componenti di STARBoard / STARBoard Komponentencodes / STARBoard-onderdeelcodes

## Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità tabella immagini / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

### Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Accertarsi che la geometria della tabella rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

<b>Table Width:</b> / Anchura de la mesa / Largeur de la table / Larghezza tavolo / Tischbreite / Tafelbreedte	45 cm - 70 cm
<b>Table Length:</b> / Longitud de la mesa / Longueur de la table / Lunghezza tavolo / Tischlänge / Tafellengte	≥200 cm
<b>Table Angulation:</b> / Angulación de la mesa / Angulation de la table / Inclinazione tavolo / Winkelstellung des Tisches / Tafelhoek	+/- 15° lateral & longitudinal +/- 15° lateral y longitudinal +/- 15° latéral et longitudinal +/- 15° laterale e longitudinale +/- 15° lateral und longitudinal +/- 15° lateraal en longitudinaal
<b>Table Surface:</b> / Superficie de la mesa / Surface de la table / Superficie del tavolo / Tischoberfläche / Tafeloppervlak	Flat Plana Plane Piatto Flach Vlak

## ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Grundlegende Informationen / Essentiële informatie

### 11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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## Information

### Intended Purpose

**STARBoard:** To support a supine patient's arm, facilitating hyperextension for radial access.

**Extension Tray:** To provide a transitional work surface from the radial access site to trolley.

**IR Shield:** To provide scatter radiation protection.

### Intended User/Training Requirement

Intended to be used by trained medical professionals.

### Patient Target Group

Adult male or female patient requiring diagnostic imaging and image guided therapy.

### Contraindications

- Not to be used with patients over 135Kg.
- Not for paediatric use.

### Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Do not use with patients over 135Kg.
- To provide restraint for restless patient (refer to page 6).
- Not for paediatric use.

### Incident Reporting Guidelines

For product complaints and incidents, please complete form on

[adeptmedical.co.nz/repairs](https://adeptmedical.co.nz/repairs)

## Initial Checks & Cleaning

### Cleaning

Refer to Disinfection Instructions.

### Product Checks

- No exposed Carbon Fibre.
- No break in Carbon Fibre surface.
- No visible damage or sharp edges, e.g. cracks.

#### STARBoard

- Ensure all pivoting joints can rotate with the appropriate amount of friction and adjust if necessary. (refer to page 5).
- Correct placement (Refer to page 5).
- Patient can access the handle to grip (refer to page 6).

#### IR Shield

- Ensure correct daggerboard placement (refer to page 7)
- Ensure the X-Ray Shield is correctly placed in the daggerboard (refer to page 7)

#### Arm Support Pad

- No visible damage or tearing.

## Disinfection

### Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.
- Some disinfectants may cause slight discolouration to the soft blue material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

### Limitations on Processing

Discontinue use if:

- Any cracks, breakages, or tears are present.
- Any exposed carbon fibre is visible.

## Disinfection Instructions

### Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

### Preparation before cleaning

**Extension Tray:** Detach the Extension Tray from the STARBoard.

**STARBoard:** Remove STARBoard from beneath the mattress.

**Arm Support Pad:** Remove Arm Support Pad from STARBoard.

**IR Shield:** Remove from Daggerboard.

**Daggerboard:** Remove from beneath the mattress.

[adeptmedicaltraining.com/productcleaning/IR-System](https://adeptmedicaltraining.com/productcleaning/IR-System)

### Cleaning: Manual

- If any visible residue is present, rinse under water, but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

### Disinfection

Refer to the Approved Cleaners List (following page) for selection of an appropriate disinfecting agent.

#### STARBoard

With a cloth dampened by an approved disinfectant, wipe all surfaces, Armrest, Handle, Wrist Pivot, and adjust device to wipe entire Slider Arm.

#### Extension Tray

Wipe all surfaces and the bracket area with a cloth dampened by an approved disinfectant.

#### IR Shield

With a cloth dampened by an approved disinfectant, wipe all surfaces including inside the slot on the Daggerboard.

#### Arm Support Pad

With a cloth dampened by an approved disinfectant, wipe all surfaces.

## Approved Cleaners List

### Approved Disinfectant by Brands

**Liquids:**

- ORION Laboratories
- 70% Isopropyl Alcohol
  
- Jaychem Industries
- 2% Chlorhexidine / 70% Alcohol
  
- Betadine
- Povidone Iodine (7.5%)
  
- BODE Chemie GmbH
- Kohrsolin FF
  
- STERIS
- Coverage Spray TB
  
- CaviCide
- Metrex Research

**Dissolvable:**

- Du Pont
- Rely+On Virkon Tablets

**Surface Wipes:**

- PDI Sani-Cloth Bleach
  
- PDI Sani-Cloth AF3
  
- Reynard Health Supplies
- Surface Disinfectant Wipes
  
- Metrex Research
- CaviWipes
  
- Clinell Wipes - Universal (Green)
  
- Clorox Hydrogen Peroxide Wipes

### Approved Active Disinfectant Ingredients

**Acid Based:**

- ≤10% Malic Acid CAS 6915-15-7
- ≤6% Sulfamic acid CAS 5329-14-6

**Alcohol Based:**

- ≤5% 2-Butoxyethanol CAS 111-76-2
- ≤10% Butyldiglycol CAS 112-34-5
- ≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0
- ≤10% Tridecanol CAS 69011-36-5
- <10% Alcohols, C12-14, ethoxylated CAS 68439-50-9
- ≤70% Denatured Ethanol CAS 64-17-5

**Sulphate Based:**

- ≤0.1% PHMB CAS 27083-27-8
- ≤55% Potassium Peroxymonosulfate CAS 70693-62-8
- ≤3% Potassium Persulfate CAS 7727-21-1

**Ammonium and Chloride Based:**

- ≤5% Benzalkonium Chloride CAS 68424-85-1
- ≤0.28% Benzethonium Chloride CAS 121-54-0
- ≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5
- ≤2% Chlorhexidine CAS 55-56-1
- ≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5
- ≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6
- <5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

**Others:**

- ≤10% Glutaral CAS 111-30-8
- ≤0.63% Sodium Hypochlorite CAS 7681-52-9
- ≤7.5% Povidone Iodine CAS 25655-41-8
- 0.5% Disodium Cocoampho Dipropionate CAS68604-71-7
- ≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via [adeptmedicaltraining.com/downloads](https://adeptmedicaltraining.com/downloads)

## Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on [adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz) for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

[adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz)

### Drying

All components should be dried thoroughly before use.

### Maintenance, Inspection and Testing

**All components:**

- Ensure no visible damage or sharp edges, e.g. cracks.

**STARBoard, X-Ray Shield:**

- Ensure no exposed Carbon Fibre is visible.

**STARBoard:**

- Ensure all pivoting joints can rotate with the appropriate amount of friction and adjust if necessary.

**Extension Tray:**

- Ensure the Extension Tray is correctly attached and secure on the STARBoard.

### Packaging

It is not required to package the device following disinfection.

### Sterilisation

This device should not be subjected to sterilisation processes.

### Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment.

*The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.*

### Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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