# PRONE SUPPORT

## AM0700

## INSTRUCTIONS FOR USE

(Instrucciones de uso / Mode d'emploi / Istruzioni per l'uso / Gebrauchsanweisung / Gebruiksaanwijzing)

Issue Date: 07/2022

Revision: B

EC REP

### EC Representative

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### Manufacturer

Adept Medical Ltd  
2-6 McDonald St,  
Morningside  
Auckland 1025, New Zealand

# CONTENTS

Contenido / Contenu / Sommario / Inhalt / Inhoud

- 4** **Unboxing** (Desembalaje / Déballage / Estrazione dalla confezione / Entpacken / Inhoud van de doos)
- 5** **Assembly** (Montaje / Assemblage / Assemblaggio / Montage / Montage)
- 6** **Setup** (Instalación/ Installation / Configurazione / Einrichtung / Installatie)
- 8** **Spare Parts** (Piezas / Pièces de rechange / Parti di ricambio / Ersatzteile / Reserveonderdelen)
- 9** **Image Table Compatibility** (Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità del tavolo per imaging / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

## **11** **Essential Information (EN)**

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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


|                                                                                    |          |                     |
|------------------------------------------------------------------------------------|----------|---------------------|
|   | <b>A</b> | Face Pad            |
|   | <b>B</b> | Face Pad Support    |
|   | <b>C</b> | Adjustable Assembly |
|  | <b>D</b> | Base Board          |

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




|                                                                                     |          |                     |       |
|-------------------------------------------------------------------------------------|----------|---------------------|-------|
|    | <b>A</b> | Face Pad            | M6531 |
|  | <b>B</b> | Face Pad Support    | M6532 |
|  | <b>C</b> | Adjustable Assembly | M6533 |
|  | <b>D</b> | Base Board          | M6534 |

### Adjustable Assembly (M6533)

## Image Table Compatibility

(Compatibilidad de la mesa de imagen / Compatibilité de la table d'examen / Compatibilità del tavolo per imaging / Kompatibilität Bildgebungstische / Compatibiliteit van MRI-tafel)

### Ensure table geometry falls within this specification

- / Asegúrese de que la geometría de la mesa cumpla estas especificaciones
- / Assurez-vous que la géométrie de la table correspond à cette spécification
- / Accertarsi che la geometria del tavolo rientri in questa specifica
- / Stellen Sie sicher, dass die Geometrie des Bildgebungstischs innerhalb dieser Angaben liegt
- / De geometrie van de tafel moet binnen deze specificatie vallen

|                                                                                                                                                         |                                                                                                                                                                                                                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Table Width:</b><br>/ Anchura de la mesa<br>/ Largeur de la table<br>/ Larghezza tavolo<br>/ Tischbreite<br>/ Tafelbreedte                           | 45 cm - 60 cm                                                                                                                                                                                                    |
| <b>Table Length:</b><br>/ Longitud de la mesa<br>/ Longueur de la table<br>/ Lunghezza tavolo<br>/ Tischlänge<br>/ Tafellengte                          | ≥200 cm                                                                                                                                                                                                          |
| <b>Table Angulation:</b><br>/ Angulación de la mesa<br>/ Angulation de la table<br>/ Inclinazione tavolo<br>/ Winkelstellung des Tisches<br>/ Tafelhoek | +/- 15° lateral & longitudinal<br>+/- 15° lateral y longitudinal<br>+/- 15° latéral et longitudinal<br>+/- 15° laterale e longitudinale<br>+/- 15° lateral und longitudinal<br>+/- 15° lateraal en longitudinaal |
| <b>Table Surface:</b><br>/ Superficie de la mesa<br>/ Surface de la table<br>/ Superficie del tavolo<br>/ Tischoberfläche<br>/ Tafeloppervlak           | Flat / Curved<br>Plana / Curva<br>Plane / incurvée<br>Piatto / Curvo<br>Flach / Gebogen<br>Vlak / gebogen                                                                                                        |
| <b>CT Bore:</b><br>/ Diámetro del túnel de TC<br>/ Alésage de l'appareil de TDM<br>/ Diametro TC<br>/ Tunneldurchmesser CT<br>/ CT-opening              | ≥78 cm                                                                                                                                                                                                           |

# ESSENTIAL INFORMATION

Información básica / Informations essentielles / Informazioni essenziali / Grundlegende Informationen / Essentiële informatie

## 11 Essential Information (EN)

- Information
- Initial Checks & Cleaning
- Disinfection
- Disinfection Instructions
- Approved Cleaners List
- Disposal

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## Information

### ? Intended Purpose

To provide a support surface for patient's arm and head in prone position.

### Intended User/Training Requirement

Intended to be used by trained medical professionals.

### Patient Target Group

Adult male or female patient undergoing procedures in prone position.

### ⚠️ Contraindications

- Not to be used with patients over 135Kg.
- Not for paediatric use.
- Not intended for CT Bore diameter less than 78cm.
- Not to be used on patients with extreme Kyphosis.
- Head Assembly not to be used with intubated patient

### ⚠️ Warnings and Cautions

- Ensure IFU is read prior to use.
- Ensure Preparatory Cleaning is conducted.
- Ensure Preparatory Product Check is conducted.
- Do not use with patients over 135Kg.
- Ensure patient has adequate mobility for shoulder flexion.
- Not for paediatric use.
- Ensure patients fingers, if over edge, do not impact with Xray camera C-arm.
- Ensure hair is managed in hair net for patients with long hair.
- If comfortable position cannot be achieved or patient is restless, the head assembly can be removed and replaced with a pillow.
- Warning, a non-cooperative patient may dislodge device as it is not fixed to Image table.
- Use of the Product is intended only for contact with intact skin.
- For patients suspected of requiring an emergency transition to the supine position during a procedure, do not use Head Assembly

### 🔧 Incident Reporting Guidelines

For product complaints and incidents, please complete form on:

[adeptmedical.co.nz/repairs](https://adeptmedical.co.nz/repairs)

## Initial Checks & Cleaning

### 🧼 Cleaning

Refer to Disinfection Instructions.

### 📋 Product Checks

- No break in Carbon Fibre surface.
- No visible damage or sharp edges, e.g. cracks
- **Adjustable Assembly:** Ensure the assembly is not free to move up and down when lever is in the lock position.
- **Adjustable Assembly:** Ensure the assembly can move up and down when pushed when lever in the unlocked position.
- **Adjustable Assembly:** Ensure the assembly is not free to tilt when the screw is locked.
- **Adjustable Assembly:** Ensure the assembly is free to tilt when the screw is unlocked.
- Ensure **Face Pad** pin is inserted securely into locating holes.
- Ensure the **Adjustable Assembly** is fully inserted on the **Face Support** (refer to page 5).
- Ensure the **Adjustable Assembly** is fully inserted on the **Base Board** (refer to page 6).

## Disinfection

### ⚠️ Warnings

- Insufficient cleaning may compromise disinfection process and lead to transmission of pathogens.
- Do not use any disinfectant products or abrasive/corrosive agents which are not on the approved cleaners list.
- Always read manufacturer's instructions and consult the manufacturers MSDS for cleaning and disinfectant products.
- Do not use any cleaning procedure which exceeds 65°C/149°F.
- Some disinfectants may cause slight discolouration to the soft blue or white material used on some components within the product range. This will not affect the strength and the product will remain fit for purpose.

### 🚫 Limitations on Processing

Discontinue use if:

- any cracks or breakages are present.
- any exposed carbon is visible.

## Disinfection Instructions

### 📋 Initial treatment at the point of use

It is important to clean the product once it is removed from the packaging and after each procedure. Ensure all areas including joints, clips, sockets, brackets and levers are thoroughly cleaned after use to remove all contaminant build-up that may be present following a procedure.

### 📋 Preparation before cleaning

Remove **Prone Support** from beneath the mattress.

Remove **Adjustable Assembly**.

Remove **Face Pad**.

Ensure the **Adjustable Assembly** is fully extended.

[adeptmedicaltraining.com/productcleaning/Prone-Support](https://adeptmedicaltraining.com/productcleaning/Prone-Support)

### 🧼 Cleaning: Manual

- If any visible residue is present, rinse under water but avoid submerging the device.
- A soft bristle brush can be used.
- Ensure the device is completely dried prior to disinfection.

### 🧴 Disinfection

Refer to the Approved Cleaners List (following page) for selection of an appropriate disinfecting agent.

#### Face Pad

With a cloth dampened by an approved disinfectant, wipe all surface.

#### Adjustable Assembly

With a cloth dampened by an approved disinfectant, wipe all surface, especially the blue lever and screw.

#### Carbon Fiber Base Board

With a cloth dampened by an approved disinfectant, wipe the entire surface, including crevices in which contaminants may build up.

### 💧 Drying

All components should be dried thoroughly before use.

### 📋 Maintenance, Inspection and Testing

For all components ensure:

- No exposed Carbon Fibre is visible.
- No visible damage or sharp edges, e.g. cracks.

#### Face Pad

- Ensure there is no damage, e.g. tear or pin detachment.

#### Adjustable Assembly

- Ensure the assembly is free to move up and down when lever is in unlocked and does not move up and down when lever is locked.

#### Adjustable Assembly

- Ensure the assembly is not free to tilt when the screw is locked and is free to tilt when the screw is unlocked.

### 📦 Packaging

It is not required to package the device following disinfection.

### 🧴 Sterilisation

This device should not be subjected to sterilisation processes.

### 📦 Storage

Once disinfection is complete and all components are dry the device should be stored in a dry environment in which the temperature will not exceed 65°C/149°F.

⚠️ *The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.*

## | Approved Cleaners List

### Approved Disinfectant by Brands

#### Liquids:

ORION Laboratories  
70% Isopropyl Alcohol

Jaychem Industries  
2% Chlorhexidine / 70% Alcohol

Betadine  
Povidone Iodine (7.5%)

BODE Chemie GmbH  
Kohrsolin FF

STERIS  
Coverage Spray TB

CaviCide  
Metrex Research

#### Dissolvable:

Du Pont  
Rely+On Virkon Tablets

#### Surface Wipes:

PDI Sani-Cloth Bleach

PDI Sani-Cloth AF3

Reynard Health Supplies  
Surface Disinfectant Wipes

Metrex Research  
CaviWipes

Clinell Wipes - Universal (Green)

Clorox Hydrogen Peroxide Wipes

### Approved Active Disinfectant Ingredients

#### Acid Based:

≤10% Malic Acid CAS 6915-15-7  
≤6% Sulfamic acid CAS 5329-14-6

#### Alcohol Based:

≤5% 2-Butoxyethanol CAS 111-76-2  
≤10% Butyldiglycol CAS 112-34-5  
≤70% Isopropyl Alcohol (Propan-2-ol) CAS 67-63-0  
≤10% Tridecanol CAS 69011-36-5  
<10% Alcohols, C12-14, ethoxylated CAS 68439-50-9  
≤70% Denatured Ethanol CAS 64-17-5

#### Sulphate Based:

≤0.1% PHMB CAS 27083-27-8  
≤55% Potassium Peroxymonosulfate CAS 70693-62-8  
≤3% Potassium Persulfate CAS 7727-21-1

#### Ammonium and Chloride Based:

≤5% Benzalkonium Chloride CAS 68424-85-1  
≤0.28% Benzethonium Chloride CAS 121-54-0  
≤10% Benzyl-C23-18-Alkyl-dimethyl Ammonium Chloride CAS 8001-54-5  
≤2% Chlorhexidine CAS 55-56-1  
≤10% Didecyl Dimethyl Ammonium chloride CAS 7173-51-5  
≤0.5% Quaternary Ammonium Compounds CAS 68956-79-6  
<5% Benzyl-C 12-18 alkyldimethylammonium chlorides CAS 63891-01-5

#### Others:

≤10% Glutaral CAS 111-30-8  
≤0.63% Sodium Hypochlorite CAS 7681-52-9  
≤7.5% Povidone Iodine CAS 25655-41-8  
0.5% Disodium Cocoampho Dipropionate CAS68604-71-7  
≤1.4% Hydrogen Peroxide CAS 7722-84-1

Access the latest Approved Cleaners list via [adeptmedicaltraining.com/downloads](http://adeptmedicaltraining.com/downloads)

## | Disposal

Used product is a biohazard, decontaminate according to instruction provided in the Disinfection Instructions and reach out to the manufacturer on [adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz) for material specification, if required. Dispose according to internal clinic policy taking into consideration local regulations.

[adeptmedical@adept.co.nz](mailto:adeptmedical@adept.co.nz)

### Serious Incidents

Any serious incident which occurs in relation to the device should be reported immediately to:

- The Competent Authority of the applicable Member State.
- The Manufacturer using the Serious Incident Reporting Form found on the Adept Medical Website.

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