

For
Economy And
Performance,

Promedica™

Auxiliary

Carts For

Video Surgical

Procedures Are

Models Of

Efficiency.



The Promedica Model 110



The Promedica Model 55

Economy

Models Loaded

With The

Features

Professionals

Push For.

Whatever video surgical procedures you're performing, wasted space is a luxury you simply can't afford. Which is why, when you need a place for a second monitor or additional equipment, Promedica auxiliary carts are the best value on wheels.


We designed them with the features that professionals tell us really count.

Like a streamlined, lighter-

weight design with the strength and durability of much larger carts. Easy maneuverability. Shelf space for supplies or additional equipment, and retainer bars to keep everything in place. All for a price that offers you the best value available in specialty carts.

The Promedica 110 Model and 55 Model auxiliary carts for video surgical procedures.

They're the models of efficiency you've been pushing for.



Model 110 Specifications

- Size: 25½" w x 23" d x 56" h
- Top Shelf: 23" w x 20¼" d (Supports 20" monitor)
- Bottom Shelf: 25½" w x 23" d
- Four Shelves: 20¼" w x 16" d
- Shelf Spacing:
 - Top: 9½"
 - Middle three: 8¼"
 - Bottom: 9½"
- Weight: 155 lbs.
- Color: Blue-Grey, Mauve
- Electrical: 6-outlet, UL approved, 15' cord on wrap bracket
- Wheels: 4" rubber castors; front two lock
- 40 lb. Ballast pan mounted on bottom of cart
- CO₂ "E" tank bracket

Order No. Model 110

Model 55 Specifications

- Size: 25" w x 23" d x 40" h
- Stainless Steel Top: 23" w x 20¼" d (Supports 20" monitor)
- Two Shelves: 20¼" w x 16" d
- Bottom Shelf: 25" w x 23" d
- Shelf Spacing:
 - Top: 12"
 - Middle: 8"
 - Bottom: 10"
- Weight: 90 lbs.
- Color: Blue-Grey, Mauve
- Wheels: 4" rubber castors; front two lock
- Electrical: 6-outlet, UL approved, 15' cord

Order No. Model 55

Manufactured By:

Promedica, Inc.
114 Douglas Road East
Oldsmar, Florida 34677
800.899.5278

Distributed By:

Medical Resources Group of North America, LLC
2312 West Altorfer Drive - Peoria, IL 61615-1810
Phone: 309-691-0534 - Fax: 309-692-3298