


Blood Exsanguination

Roll-On Cuffs

- Exsanguinator allows fast and easy blood exsanguination within seconds (as alternative to Esmarch Bandage)
- Material: Silicone, transparent
- Reusable, autoclavable up to 134°C (273°F)
- Latex free

Size	Color code	Patient limb		REF
0	grey	16-21cm	6"-8"	21-95-710
1	white	20-28cm	8"-11"	21-95-711
2	light-blue	26-33cm	10"-13"	21-95-712
3	yellow	31-40cm	12"-16"	21-95-719
4	red	38-49cm	15"-19"	21-95-722
5	green	47-60cm	18"-23"	21-95-727
6	dark-blue	58-70cm	23"-28"	21-95-729



Handpump
with manometer and
1m (40") connecting tube

REF 21-95-350




Manometer
with 1m (40") connecting tube

REF 21-95-300

Brake-Block

made of Silicone,
autoclavable up to 134°C (273°F)

REF 21-95-200 small size
REF 21-95-201 large size



Tape-Measure

color coded, for easy identification
of patient extremity and cuff size

REF 21-95-111



REF 21-95-400

Roll-On Cuff Set

small sizes

complete case with Roll-On Cuff sizes 1-4, manometer with connecting tube, tape-measure, brake-block small and large

REF 21-95-400

Roll-On Cuff Set

large sizes

complete case with Roll-On Cuff sizes 3-6, handpump with manometer and connecting tube, tape-measure, brake-block small and large

REF 21-95-401

Esmarch Bandages

autoclavable up to 134°C (273°F)
 Length: 5 m (16ft.)
 latex free

Width	inch	REF
6cm	2 ³ / ₈ "	20-20-060
8cm	3 ¹ / ₈ "	20-20-080
10cm	4"	20-20-100
12cm	4 ³ / ₄ "	20-20-120
15cm	6"	20-20-150



disposable, sterile packed, Box of 10
 Length: 4 m (13ft.)
 latex free

Width	inch	REF
8cm	3 ¹ / ₈ "	20-21-080-1
10cm	4"	20-21-100-1
15cm	6"	20-21-150-1



Three Chamber Bloodexsanguination Garment

Automatic exsanguination for use with VBM Tourniquets 2800, 5000 and 5800

Due to the Three Chamber-System the garment is inflated automatically from distal to proximal.

With the 3 position zipper the garment can be adjusted to fit different size extremities for efficient exsanguination of the patients limb. Sleeve and bladder are exchangeable and available as spare parts.



REF 20-13-341



REF 20-13-343

Three Chamber Bloodexsanguination Garment, for arm
 Three Chamber Bloodexsanguination Garment, for leg

REF 20-13-341
 REF 20-13-343

