

# Tourniquet Cuffs

## Product details

### Silicone Bladder

- Silicone material is elastic and acts as a low pressure cuff (as wider the cuff as lower the inflation pressures needed)
- Durable and high quality material
- Latex free

### Cuff ties

- Provide easy cuff positioning and security
- Color coded for identification of different sizes

### Silicone Tubing

- Kink resistant
- Never cracks
- Latex free

### Positive locking connectors (PLC)

- For easy attachment and detachment
- For safe and leak-free inflation

## Single Cuffs



### Fabric Cuffs

- Washable and reusable
- Fixation via Velcro tape
- Bladder and fabric sleeve can be replaced



### Silicone Cuffs

- Autoclavable at 134°C (273°F)
- Fixation via strap and buckle
- Easy to clean and dry quickly



### Silicone/Velcro Cuffs

- Autoclavable at 134°C (273°F)
- Fixation via Velcro tape
- Velcro sleeve can be replaced
- Easy to clean and dry quickly



### Disposable Cuffs

- Disposable, sterile peel pouch
- Low profile, welded cuff (no bladder)
- Eliminate risk of cross contamination


Patient,  
Size

				PLC Connector	Luer Lock Connector
Infant, 20cm (8")	20-54-700	—	20-64-700	20-34-700SLZ-1	20-34-700LL-1
Child, 30cm (12")	20-54-710	20-94-710	20-64-710	20-34-710SLZ-1	20-34-710LL-1
Arm, 35cm (14")	20-54-711	20-94-711	20-64-711	20-34-711SLZ-1	20-34-711LL-1
Arm, 46cm (18")	20-54-712	20-94-712	20-64-712	20-34-712SLZ-1	20-34-712LL-1
Leg, 61cm (24")	20-54-722	20-94-722	20-64-722	20-34-722SLZ-1	20-34-722LL-1
Leg, 76cm (30")	20-54-727	20-94-727	20-64-727	20-34-727SLZ-1	20-34-727LL-1
Leg, 86cm (34")	20-54-728	20-94-728	20-64-728	20-34-728SLZ-1	20-34-728LL-1
Leg, 107cm (42")	20-54-729	20-94-729	20-64-729	20-34-729SLZ-1	20-34-729LL-1

## Plastic Cover

### Tourniquet Cuff Protector

The plastic covers keep the cuff free from soiling and allow sealing the cuff with adhesive tape without damaging the cuff.



REF	Size	Box
REF 20-11-111	small	50
REF 20-11-111-1	ditto but sterile	50
REF 20-11-222	large	50
REF 20-11-222-1	ditto but sterile	50

## Tape-Measure

color coded, for easy identification of patient extremity and cuff size

REF 20-95-111



## Double Cuffs



### Fabric Cuffs

- Washable and reusable
- Fixation via Velcro tape
- Bladder and fabric sleeve can be replaced



### Silicone Cuffs

- Autoclavable at 134°C (273°F)
- Fixation via strap and buckle
- Easy to clean and dry quickly



### Silicone/Velcro Cuffs

- Autoclavable at 134°C (273°F)
- Fixation via Velcro tape
- Velcro sleeve can be replaced
- Easy to clean and dry quickly

Patient, Size	Fabric Cuffs	Silicone Cuffs	Silicone/Velcro Cuffs
Infant, 20cm (8")	20-50-700	—	20-60-700
Child, 30cm (12")	20-50-710	20-90-710	20-60-710
Arm, 35cm (14")	20-50-711	20-90-711	20-60-711
Arm, 46cm (18")	20-50-712	20-90-712	20-60-712
Leg, 61cm (24")	20-50-722	20-90-722	20-60-722
Leg, 76cm (30")	20-50-727	—	20-60-727
Leg, 86cm (34")	20-50-728	—	—
Leg, 107cm (42")	20-50-729	—	—

