

# ENT OFFICE-BASED RADIO-FREQUENCY SURGICAL PROCEDURES




## ESU 70B SURGICAL PROCEDURE SYSTEM

MADE IN THE USA

When safety, simplicity, service-free operation, and economical pricing are your objective, the **ELMED ESU 70 B** and **ESU 120 M/M DIGITAL** are the units for you. **ELMED**, the leader in electrosurgery manufactures its units and reusable accessories completely in-house, in Glendale Heights, IL, USA with the aim to reduce your surgical cost.

Both units and reusable accessories are specifically designed for all office-based otolaryngology procedures, including monopolar and bipolar surgical procedures.



The ELMED ESU 70 B Provides Treatment for:

- Nose Bleeds
  - Treatment of Snoring
  - Submucosal Turbinate Surgery
  - Obstructive Tonsillar Hypertrophy
- All monopolar/bipolar cut, biopsy & coagulation procedures

LIST NO.	ESU 70 B SURGICAL PROCEDURE SYSTEM COSISTING OF ALL THE ITEMS LISTED BELOW
----------	--

5196	ESU 70 B (includes footswitch)
5047*	PATIENT PLATE
5701*	PATIENT PLATE CABLE
5044*	RUBBER STRAP
5021	MONOPOLAR RF-CABLE
5241-1	SUCTION COAG. HANDLE & EXTENSION
5004	FOOTSWITCH HANDLE
5807-05	BLADE ELECTRODE
5970	DENNIS PROBE
5230	BIPOLAR RF-CABLE
4755	EQUIPMENT TROLLEY
5798B-1	RFTA HANDLE & ELECTRODE

\*Disposable patient plates available on request



## ESU 120 M/M DIGITAL SURGICAL PROCEDURE SYSTEM

MADE IN THE USA



The ELMED ESU 120 M/M Digital Provides Treatment for:

- Nose Bleeds
  - Treatment of Snoring
  - Submucosal Turbinate Surgery
  - Obstructive Tonsillar Hypertrophy
- All monopolar/bipolar cut, biopsy & coagulation procedures

LIST NO.	ESU 120 M/M DIGITAL SURGICAL PROCEDURE SYSTEM COSISTING OF ALL THE ITEMS LISTED BELOW
5797 ENT II	

5797	ESU 70 B (includes footswitch)
5047*	PATIENT PLATE
5701*	PATIENT PLATE CABLE
5044*	RUBBER STRAP
5021	MONOPOLAR RF-CABLE
5241-1	SUCTION COAG. HANDLE & EXTENSION
5004	FOOTSWITCH HANDLE
5807-05	BLADE ELECTRODE
5970	DENNIS PROBE
5302	BIPOLAR RF-CABLE
4755	EQUIPMENT TROLLEY
5798B-1	RFTA HANDLE & ELECTRODE

\*Disposable patient plates available on request