


HEGA - MEDICAL

Cleaning Technology for Medical and Research Facilities



1925 South Main Street • Morton, IL 61550 • O: 309.263.6901 • T: 888.674.6621
F: 309.263.6903 • E: sales@mrgona.com • www.mrgona.com

MEDICAL RESOURCES GROUP LIMITED

Developed from experience, **HEGA** Rinsing Guns stand for clean results with safe and easy handling.

Whether for cleaning instrumentation, endoscopes, cystoscopies, laboratory glassware, or any lumened instrument, **HEGA** Rinsing and Washing Guns are the perfect choice.


HEGA Rinsing Guns ensure high volume regardless of drying with air or cleaning with water. The **HEGA** accessory tips ensure every remote corner is reached.

- Up to 145psi
- Use with Deionized Water
- Chrome Plated Brass for Corrosion Resistance
- Attachment Holder Polished Stainless Steel





TERMS and Conditions:



- 1.) Terms are Net 30 days
- 2.) Shipping is FOB Ship Point Pre Pay and Add
- 3.) All HEGA-Medical Products are guaranteed for material and workmanship for the period from 12 months from date of purchase.
- 4.) Returns will be accepted 60 days from date of purchase.
- 5.) We cannot accept the following for credit or refund:
 - a. Instruments not in their original condition
 - b. Instruments that have been used, unless defective. Used instruments must be decontaminated and include documentation of decontamination.
- 6.) Prices are subject to change without notice.



Dina Rinsing Gun with Mountable Connection Set



Dina Rinsing Gun with Mountable Connection Set				
Product Number	Description	Connection	UOM	
6015	1.5M (59.06") PVC Fabric Reinforced Hose	1/2"	EA/1	
6025	1.5M (59.06") PVC Fabric Reinforced Hose	3/4"	EA/1	
6035	3.0M (118.11") PVC Fabric Reinforced Hose	1/2"	EA/1	
6045	3.0M (118.11") PVC Fabric Reinforced Hose	3/4"	EA/1	
Product Number	Description	Connection	UOM	
6055	1.5M (59.06") Stainless Steel Hose	1/2"	EA/1	
6065	1.5M (59.06") Stainless Steel Hose	3/4"	EA/1	
6075	3.0M (118.11") Stainless Steel Hose	1/2"	EA/1	
6085	3.0M (118.11") Stainless Steel Hose	3/4"	EA/1	

Nova Rinsing Gun with Mountable Connection Set

Nova Rinsing Gun with Mountable Connection Set				
Product Number	Description	Connection	UOM	
7015	1.5M (59.06") PVC Fabric Reinforced Hose	1/2"	EA/1	
7025	1.5M (59.06") PVC Fabric Reinforced Hose	3/4"	EA/1	
7035	3.0M (118.11") PVC Fabric Reinforced Hose	1/2"	EA/1	
7045	3.0M (118.11") PVC Fabric Reinforced Hose	3/4"	EA/1	
Product Number	Description	Connection	UOM	
7055	1.5M (59.06") Stainless Steel Hose	1/2"	EA/1	
7065	1.5M (59.06") Stainless Steel Hose	3/4"	EA/1	
7075	3.0M (118.11") Stainless Steel Hose	1/2"	EA/1	
7085	3.0M (118.11") Stainless Steel Hose	3/4"	EA/1	



Dina Cleaning Kit

Dina Rinsing Gun with Professional Set Including 8 Attachments and Holder

Dina Rinsing Gun with Professional Set Including 8 Attachments and Holder				
Product Number	Description	Connection	UOM	
6010	1.5M (59.06") PVC Fabric Reinforced Hose	1/2"	EA/1	
6020	1.5M (59.06") PVC Fabric Reinforced Hose	3/4"	EA/1	
6030	3.0M (118.11") PVC Fabric Reinforced Hose	1/2"	EA/1	
6040	3.0M (118.11") PVC Fabric Reinforced Hose	3/4"	EA/1	
Product Number	Description	Connection	UOM	
6050	1.5M (59.06") Stainless Steel Hose	1/2"	EA/1	
6060	1.5M (59.06") Stainless Steel Hose	3/4"	EA/1	
6070	3.0M (118.11") Stainless Steel Hose	1/2"	EA/1	
6080	3.0M (118.11") Stainless Steel Hose	3/4"	EA/1	


Nova Cleaning Kit

Nova Rinsing Gun Professional Set Including 8 Attachments and Attachment Holder

Nova Rinsing Gun Professional Set Including 8 Attachments and Holder				
Product Number	Description	Connection	UOM	
7010	1.5M (59.06") PVC Fabric Reinforced Hose	1/2"	EA/1	
7020	1.5M (59.06") PVC Fabric Reinforced Hose	3/4"	EA/1	
7030	3.0M (118.11") PVC Fabric Reinforced Hose	1/2"	EA/1	
7040	3.0M (118.11") PVC Fabric Reinforced Hose	3/4"	EA/1	
Product Number	Description	Connection	UOM	
7050	1.5M (59.06") Stainless Steel Hose	1/2"	EA/1	
7060	1.5M (59.06") Stainless Steel Hose	3/4"	EA/1	
7070	3.0M (118.11") Stainless Steel Hose	1/2"	EA/1	
7080	3.0M (118.11") Stainless Steel Hose	3/4"	EA/1	




PG 3



1925 South Main Street • Morton, IL 61550 • O: 309.263.6901 • T: 888.674.6621
F: 309.263.6903 • E: sales@mrgona.com • www.mrgona.com

MEDICAL RESOURCES GROUP LIMITED

Rinsing Gun Only


Dina Rinsing Gun with Nozzle and Spout Only			
Product Number	Description	UOM	
6000	Dina Rinsing Gun Only	EA/1	
Nova Rinsing Gun with Nozzle and Spout Only			
Product Number	Description	UOM	
7000	Nova Rinsing Gun Only	EA/1	

Cleaning Kit Attachments

Cleaning Attachments for Dina and Nova Rinsing Guns		
Product Number	Description	UOM
4230	Universal Attachment	EA/1
4240	Lumen Tip Attachment	EA/1
4250	Spray Jet Attachment	EA/1
4260	Syringe and Cannula Luer Tip Attachment	EA/1
4270	Syringe and Cannula Tip Attachment	EA/1
5010	Large Lumen Attachment	EA/1
5020	Bottle and Flask Attachment	EA/1
5030	Water Jet / Suction Attachment	EA/1

4230	4240	4250	4260
			
4270	5010	5020	5030
			

HEGA Rinsing Gun Accessories

Dina Rinsing Gun Only		
Product Number	Description	UOM
5220	Professional Attachment Set: Includes 8 Attachments and 1 Attachment Holder. Holder Specifications: Stainless Steel 2mm x 280mm x 40mm, Self Adhesive	EA/1
		

Accessories		
Product Number	Description	UOM
5040	8 Attachment Holder, SS, 2mm x 280mm x 40mm	EA/1
3000	Wall Holder for HEGA Gun Polyamid	EA/1
3005	Wall Holder for HEGA Gun Chrome Plated Brass	EA/1
3020	Wall Holder for HEGA Gun Stainless Steel Hook	EA/1
4550	Silicone O-Ring	PK/3
5770	Table Grommet with Lock Nut, up to 50mm Thickness	EA/1

5040		3000	3005
			
3020	4550		5770
			