




SteriProtect

INSTRUMENT TRAY CORNER PROTECTORS

SteriProtect Corner Protectors are a simple and effective solution to eliminate sterilization wrap perforation and prevent Instrument Tray corner damage during sterilization.

Procedure trays of surgical instruments are often wrapped in non-woven CSR wrap to maintain sterility until use. If the wrap has any holes, becomes ripped or torn then it must be rejected in the operating room causing serious, untimely and costly delays.

- **SteriProtect Tray Corner Protectors** safeguard surgical instrument tray wrappers from rips, holes and tears due to sharp tray corners, improper handling or heavy surgical trays.
- **SteriProtect Tray Corner Protectors** are designed and manufactured in a soft touch TPE to create a flexible “bumper” between the tray’s sharp corners and the sterilization wrap.
- **SteriProtect Tray Corner Protectors** open cell design also allows for air to vent freely.
- **SteriProtect Tray Corner Protectors** can be used with steam, ethylene oxide and low temperature plasma sterilization.
- **SteriProtect Tray Corner Protectors** are reusable, economical to use, our pink color contrasts with blue wrappers ensuring return to Sterile Processing for reuse.
- **SteriProtect Tray Corner Protectors** are Latex and BPA free.



Catalog Number	Description	UOM
MRG0400CP	SteriProtect Instrument Tray Protector	CS/400

MEDICAL RESOURCES GROUP LIMITED