





Myringotomy Kits

Our all-in-one procedure packs contain sterilised, easy to use, single-use instruments that are economically packaged to remove the possibility of cross patient contamination and are priced to save on cleaning and assembly costs. All products are supplied sterile and individually pouched.

- Easy to use and set up.
- High-quality single-use instruments.
- Ready to use.
- Sterilised by ETO gas.

<p>NZ3246-5 Myringotomy Kit</p>	<p>5 / Pack</p>
<p>1 x Myringotomy Knife</p>	
<p>1 x Jobson Horne Probe</p>	
<p>1 x Rosen Inserter</p>	
<p>2 x Bent Suction Tubes (16G & 20G)</p>	
<p>1 x Suction Controller</p>	
<p>3 x Specula (4mm, 5mm, 6mm)</p>	
<p>NZ3248B-5 Myringotomy Kit - 7mm Jaw</p>	<p>5 / Pack</p>
<p>1 x Myringotomy Knife</p>	
<p>1 x Jobson Horne Probe</p>	
<p>1 x Forceps - Stainless Steel, 7mm Jaw</p>	
<p>1 x Rosen Inserter</p>	
<p>2 x Bent Suction Tubes (16G & 20G)</p>	
<p>1 x Suction Controller</p>	
<p>3 x Specula (4mm, 5mm, 6mm)</p>	
<p>NZ3270-5 Myringotomy Kit - 3.5mm Jaw</p>	<p>5 / Pack</p>
<p>1 x Myringotomy Knife</p>	
<p>1 x Jobson Horne Probe</p>	
<p>1 x Forceps - Stainless Steel, 3.5mm Jaw</p>	
<p>1 x Rosen Inserter</p>	
<p>2 x Bent Suction Tubes (16G & 20G)</p>	
<p>1 x Suction Controller</p>	
<p>3 x Specula (4mm, 5mm, 6mm)</p>	
<p>NZ3271B-5 Myringotomy Kit Paediatric</p>	<p>5 / Pack</p>
<p>1 x Myringotomy Knife</p>	
<p>1 x Jobson Horne Probe, Oval</p>	
<p>1 x Ultra Fine Forceps - Stainless Steel, 3.5mm Jaw</p>	
<p>1 x Rosen Inserter</p>	
<p>2 x Bent Suction Tubes (16G & 18G)</p>	
<p>1 x Suction Controller</p>	
<p>3 x Ear Specula (3mm, 4mm, 5mm)</p>	