

Otological Instruments

A comprehensive range of single-use myringotomy instruments for surgical theatres and clinics.

NeoZoline™ offers innovative cerumen management products for hospital, primary care or audiology clinic use. Our single-use products perform better, are safer and more cost-efficient than traditional metal instruments.

CERUMEN HOOK

Unique single-use Cerumen Hook allows ease of loosening hard cerumen (ear wax) and other obstacles from the ear canal.

- Optimised hook geometry gives maximum leverage.
- Reduced patient trauma due to smooth edges.

- Moulded from carbon fibre reinforced nylon for strength and rigidity, resulting in a light-weight, flexible and strong single-use version of the traditional metal instrument.



Code	Description	Material	Quantity
NZ3441-100	Cerumen Hook, Non-Sterile	Carbon Fibre Filled Nylon	100 / Pack

JOBSON HORNE PROBE - EAR CURETTE

Designed to allow gentle and easy extraction of loosened cerumen (ear wax) from the ear canal rather than pushing it deeper inside.

- Moulded from carbon fibre reinforced nylon for strength and rigidity.
- Optimally angled for safe and effective removal of loosened cerumen (ear wax).

- Probe is designed as a dual purpose instrument with a curette loop at one end and a threaded section at the other to secure wrapped cotton wool.



Code	Description	Loop Diameter	Material	Quantity
NZ3431-100	Jobson Horne Probe, Non-Sterile	3 mm	Carbon Fibre Filled Nylon	100 / Pack
NZ3431-20S	Jobson Horne Probe, Sterile *	3 mm	Carbon Fibre Filled Nylon	20 / Pack
NZ3432-100	Jobson Horne Probe - Oval, Non-Sterile	2.5 mm	Carbon Fibre Filled Nylon	100 / Pack
NZ3432-20S	Jobson Horne Probe - Oval, Sterile	2.5 mm	Carbon Fibre Filled Nylon	20 / Pack

ROSEN INSERTER

- Optimised hook geometry gives maximum leverage.
- Reduced patient trauma due to smooth edges.

- Moulded from carbon fibre reinforced nylon for strength and rigidity, resulting in a light-weight, flexible and strong single-use version of the traditional metal instrument.



Code	Description	Material	Quantity
NZ3461-100	Rosen Inserter, Non-Sterile	Carbon Fibre Filled Nylon	100 / Pack

MYRINGOTOMY KNIFE

NeoZoline™ Myringotomy Knife is used to pierce and incise the tympanic membrane. The blade is sharpened and angled to aid incision.

- Raised area on handle for touch orientation.
- Bendable shaft for ease of use, versatility and prevents obstruction of the optic axis.

- Made of medical grade stainless steel and polypropylene, making it compliant with ISO 10993.



Code	Description	Length	Material	Quantity
NZ3245-10	Myringotomy Knife - Lancet, Sterile *	157mm	Stainless Steel & Polypropylene	10 / Pack

Otological Instruments

MICRO EAR SUCTION TUBES

Suction Tubes for suction of fluid and other media from the ear.

- Suction Tube has colour coded hub for ease of tube size identification.
- Cost-efficient and contamination risk-free alternative to reusable versions.
- Easy to use.
- Bent to an optimal angle.



Code	Description	OD	Colour	Quantity
NZ4254Bent	Suction Tube 14 Gauge Bent, Sterile *	2.1 mm	Green	20 / Pack
NZ4250Bent	Suction Tube 16 Gauge Bent, Sterile *	1.6 mm	White	20 / Pack
NZ4251Bent	Suction Tube 18 Gauge Bent, Sterile *	1.2 mm	Pink	20 / Pack
NZ4252Bent	Suction Tube 20 Gauge Bent, Sterile *	0.9 mm	Yellow	20 / Pack
NZ4253Bent	Suction Tube 22 Gauge Bent, Sterile *	0.7 mm	Grey	20 / Pack
NZ4255Bent	Suction Tube 24 Gauge Bent, Sterile *	0.5 mm	Purple	20 / Pack
NZ5254Bent	Suction Tube 14 Gauge Bent 100 Pack, Non-Sterile	2.1 mm	Green	100 / Pack
NZ5250Bent	Suction Tube 16 Gauge Bent 100 Pack, Non-Sterile	1.6 mm	White	100 / Pack
NZ5251Bent	Suction Tube 18 Gauge Bent 100 Pack, Non-Sterile	1.2 mm	Pink	100 / Pack
NZ5252Bent	Suction Tube 20 Gauge Bent 100 Pack, Non-Sterile	0.9 mm	Yellow	100 / Pack
NZ5253Bent	Suction Tube 22 Gauge Bent 100 Pack, Non-Sterile	0.7 mm	Grey	100 / Pack
NZ5255Bent	Suction Tube 24 Gauge Bent 100 Pack, Non-Sterile	0.5 mm	Purple	100 / Pack

MICRO EAR SUCTION TUBE AND CONTROLLER PACKS

Suction Controller governs level of suction through a thumb vent mechanism.

- Colour coded hubs for easy identification of tube size
- Easy to use.
- Easy to assemble.
- Bent to an optimal angle.
- Pack contains tube and controller.



Code	Description	OD	Colour	Quantity
NZ3254Bent	Suction Tube 14 Gauge Bent, Sterile *	2.1 mm	Green	10 / Pack
NZ3250Bent	Suction Tube 16 Gauge Bent, Sterile *	1.6 mm	White	10 / Pack
NZ3251Bent	Suction Tube 18 Gauge Bent, Sterile *	1.2 mm	Pink	10 / Pack
NZ3252Bent	Suction Tube 20 Gauge Bent, Sterile *	0.9 mm	Yellow	10 / Pack
NZ3253Bent	Suction Tube 22 Gauge Bent, Sterile *	0.7 mm	Grey	10 / Pack
NZ3255Bent	Suction Tube 24 Gauge Bent, Sterile *	0.5 mm	Purple	10 / Pack

* Individually pouched – per pack.

Otological Instruments

SUCTION CONTROLLER

The Suction Controller governs the level of suction through a thumb vent mechanism.

- Designed to fit NeoZoline™ suction tubes.
- Easy to use.



Code	Description	Colour	Length	Material	Qty
NZ3238-100	Suction Controller, Non-Sterile	Black	39 mm	ABS	100 / Pack
NZ3238S	Suction Controller, Sterile *	Black	39 mm	ABS	20 / Pack

EAR SPECULA

Advanced design provides maximum patient comfort and superior operational vision.

- Matte finish and internal ridges reduce glare and reflection on the internal surface of the speculum from the microscope light.
- Six medial end sizes are available in non-sterile polyethylene.
- Slightly oval in shape to conform to the ear canal.
- Soft and flexible material to aid patient comfort.
- Serrated lateral end to aid finger grip.



Code	Description	Material	Quantity
NZ3235-5	Ear Specula 10mm, Non-Sterile	Polyethylene	5 / Pack
NZ3234-5	Ear Specula 8mm, Non-Sterile	Polyethylene	5 / Pack
NZ3230-5	Ear Specula 6mm, Non-Sterile	Polyethylene	5 / Pack
NZ3231-5	Ear Specula 5mm, Non-Sterile	Polyethylene	5 / Pack
NZ3232-5	Ear Specula 4mm, Non-Sterile	Polyethylene	5 / Pack
NZ3233-5	Ear Specula 3mm, Non-Sterile	Polyethylene	5 / Pack
NZ3235-100	Ear Specula 10mm, Non-Sterile	Polyethylene	100 / Pack
NZ3234-100	Ear Specula 8mm, Non-Sterile	Polyethylene	100 / Pack
NZ3230-100	Ear Specula 6mm, Non-Sterile	Polyethylene	100 / Pack
NZ3231-100	Ear Specula 5mm, Non-Sterile	Polyethylene	100 / Pack
NZ3232-100	Ear Specula 4mm, Non-Sterile	Polyethylene	100 / Pack
NZ3233-100	Ear Specula 3mm, Non-Sterile	Polyethylene	100 / Pack

MICRO EAR FORCEPS

Single-use stainless steel Micro Ear Forceps.

- Made of medical grade stainless steel, making them compliant with ISO 10993.
- Tapered tip for ease of manoeuvrability in the ear canal.
- Three variations of size to suit a wide range of patients.
- Smooth opening and closing action.
- Comfortable to use.



Code	Description	Jaw Length	Length	Qty
NZ3452-10NS	Micro Ear Forceps, Stainless Steel, Non-Sterile	3.5 mm	75 mm	10 / Pack
NZ3453B-10NS	Micro Ear Forceps, Stainless Steel, Non-Sterile	7 mm	85 mm	10 / Pack
NZ3457-10NS	Ultra Fine Ear Forceps, Stainless Steel, Non-Sterile	3.5 mm	75 mm	10 / Pack