

Rhinology

SELF-ADHESIVE EXTERNAL NASAL SPLINT

Following Rhinoplasty procedures, the splint provides protection for the nose against trauma as well as preventing excessive soft tissue swelling which may distort nasal bone alignment.

- Light-weight aluminium which moulds to any shape of nose and can be trimmed easily to adjust size, but remains rigid once applied.
- Dorsal sponge support for patient comfort.
- Adhesive foam strips allow for easy and efficient application and removal.



Code	Description	Dimensions	Colour	Material	Qty
NZ3243_SML	Nasal Splint—External, Small, Non-Sterile	32mm x 56mm	Almond	Aluminium	20 / Pack
NZ3243	Nasal Splint—External, Medium, Non-Sterile	43mm x 62mm	Almond	Aluminium	20 / Pack
NZ3243_LGE	Nasal Splint—External, Large, Non-Sterile	55mm x 75mm	Almond	Aluminium	20 / Pack

TILLEY NASAL FORCEPS

- Manufactured from medical grade stainless steel.
- Cost-efficient and contamination risk-free alternative to reusable versions.



Code	Description	Material	Qty
NZ6610-10NS	Tilley Nasal Forcep, 140mm, Non-Sterile	Stainless Steel	10 / Pack

THUDICUM NASAL SPECULA

- Available in three blade sizes allowing for selection of speculum appropriate to the patient.
- Manufactured from medical grade stainless steel.



Code	Description	Material	Qty
NZ6812	Thudicum Nasal Specula Size 2, Non-Sterile	Stainless Steel	10 / Pack
NZ6814	Thudicum Nasal Specula Size 4, Non-Sterile	Stainless Steel	10 / Pack
NZ6816	Thudicum Nasal Specula Size 6, Non-Sterile	Stainless Steel	10 / Pack

Rhinology

SINUS CANNULAE WITH CONTROLLER AND NASAL SUCTION TUBES WITH CONTROLLER

Sinus Cannulae and Nasal Suction Tubes have colour coded hubs for easy identification of tube size. The Suction Controller has an easy thumb vent control mechanism and it offers a single-use, cost efficient and contamination risk-free alternative to reusable versions.



Code	Description	OD	Colour	Length	Quantity
NZ6212	Nasal Suction Tube 12 Gauge + Controller, Sterile *	2.4 mm	Red	150 mm	10 / Pack
NZ6214	Nasal Suction Tube 14 Gauge + Controller, Sterile *	2.1 mm	Blue	150 mm	10 / Pack
NZ6412	Sinus Cannulae 12 Gauge + Controller, Sterile *	2.4 mm	Red	150 mm	10 / Pack
NZ6414	Sinus Cannulae 14 Gauge + Controller, Sterile *	2.1 mm	Blue	150 mm	10 / Pack