










Otological Ventilation Tubes


Also known as Tympanostomy Tubes, Ear Tubes and Grommets, they are used in the treatment of Otitis Media. Our extensive range includes fluoroplastic and silicone tubes supplied in our easy to use dispenser.

All tubes are supplied sterile in an ergonomic case designed to ensure the tube stays secure. Are available as singles or doubles. All products are supplied sterile and individually pouched.

Code	Internal Diameter	Inner Flange Diameter	Outer Flange Diameter	Inter Flange Distance	Length	Feature	Material
Armstrong							
NZ3121	1.14	2.5	2.5	-	3.8	White	Fluoroplastic
NZ3122	1.14	2.5	2.5	-	3.8	Tail	Fluoroplastic
Modified Armstrong							
NZ3123	1.14	2.5	2.5	-	2.3	White	Fluoroplastic
NZ3124	1.14	2.5	2.5	-	2.3	Tail	Fluoroplastic
Bevel Bobbin							
NZ3141	1.14	2.8	2.8	0.7	1.7	White	Fluoroplastic
Straight Tube							
NZ3151	1.14	2.7	1.55	-	7.25	White	Fluoroplastic
Collar Button							
NZ3211	1.14	2.6	2.6	1.05	1.65	White	Fluoroplastic
NZ3211Blue	1.14	2.6	2.6	1.05	1.65	Blue	Fluoroplastic
NZ3212	1.14	2.6	2.6	1.05	1.65	Tail	Fluoroplastic
NZ3111	1.27	3.0	3.0	1.55	2.15	White	Fluoroplastic
NZ3112	1.27	3.0	3.0	1.55	2.15	Tail	Fluoroplastic

Otological Ventilation Tubes

Code	Internal Diameter	Inner Flange Diameter	Outer Flange Diameter	Inter Flange Distance	Length	Feature	Material
							
Shepard							
NZ3213	1.14	2.4	2.4	-	2.2	White	Fluoroplastic
NZ3213Blue	1.14	2.4	2.4	-	2.2	Blue	Fluoroplastic
NZ3214	1.14	2.4	2.4	-	2.2	Tail	Fluoroplastic
							
Shah and Mini Shah							
NZ3223	0.76	2.2 *	1.5	-	1.5	White	Fluoroplastic
NZ3223Blue	0.76	2.2 *	1.5	-	1.5	Blue	Fluoroplastic
NZ3226	0.76	2.2 *	1.5	-	1.5	Tail	Fluoroplastic
NZ3224	1.14	3.4 *	2.3	-	2.4	White	Fluoroplastic
NZ3227	1.14	3.4 *	2.3	-	2.4	Tail	Fluoroplastic
*Maximum external flange dimensions							

T-Tube							
							
NZ3306	1.3	8.0	1.7	-	6.0	Blue	Silicone
NZ3309	1.3	8.0	1.7	-	9.0	Blue	Silicone
NZ3312	1.3	8.0	1.7	-	12.0	Blue	Silicone

Donaldson							
							
NZ3321	1.1	2.3	2.3	1.0	2.2	Blue	Silicone

Otological Display Case

The Otological Ventilation Tube Display Case can be used for display purposes for doctors showcasing Adept Medical's NeoZoline™ Otological Ventilation Tube Range. Educating the patient or their family members on the procedure they are having helps to ease apprehension.

The display case comes with our grommets securely captured.



Code	Description	Size	Quantity
NZM001	NeoZoline grommet range	92mm x 54mm x 8mm	1 Case